



The **Ecolab ORS-2075R** solution warmer is a fluid warming system. It provides immediate access to warm irrigation fluid at a visible and controlled temperature to increase patient safety and reduce the risk of hypothermia. The Ecolab ORS-2075R temperature settings can be adjusted and locked at a range from room temperature to 120°F (49°C).

Features

- Optional Portable Stand.
- Basin Capacity: 7.5 Liters.
- Dimensions: 14.25" O.D. Round x 4".
- Minimum Fluid Level: 1.5 liters.

