



The **Endostat III** is designed for compatibility with Injection Gold Probe™ and Gold Probe™ Electrohemostasis Catheters, polypectomy snares, sphincterotomes, and hot biopsy forceps.

Features

- A self-limiting energy system that is designed to limit the depth of burn penetration.
- A variable speed peristaltic pump intended to provide precisely controlled irrigation.
- Two bipolar modes: CONTINUOUS or VARIABLE TIME, designed to vary the duration of energy delivery.
- The monopolar capabilities of the Endostat III Generator are intended to provide physicians the ability to perform certain therapeutic procedures using standard monopolar devices such as polypectomy snares, sphincterotomes, and hot biopsy forceps.
- A choice of four individual RF waveforms: COAG, BLEND, CONTROL CUT, or CUT.
- Receptacles that accept widely used monopolar active cords and disposable dispersive electrodes (patient plates, or grounding pads).

