



Video endoscopy is a process by which the pathologist will determine the swallowing capacity of a patient. The **Fujinon EG- 450WR5** is a high-resolution camera that allows for superb visualization and crystal clear dramatic images with stunning contrast with its broad view of 140 degrees. The device can be inserted without any problems with its 2.8mm working channel. It has a thin 9.44mm outer diameter.

## Features

- Ultra-wide 140° field of view
- Fujinon's traditional 410,000 pixel CCD chip
- Incredible visualization and dramatic color replication
- Slim 9.4mm outer diameter
- 2.8mm working channel
- Easy insertion
- Standard accessory passage



# Specifications

## Optical System

**Field of View:** 120°

**Direction of View:** 0° Forward viewing

**Depth of Field:** 4 to 100 mm

## Insertion Tube

**Outer Diameter:** 9.3 mm

**Working Length:** 1100 mm

**Total Length:** 1400 mm

**Instrument Channel Size:** 2.8 mm

## Bending Section

### Angulation Range:

- **Up:** 210°
- **Down:** 90°
- **Right:** 100°
- **Left:** 100°

