



The **GE 3S-RS Transducer** is compatible with the GE Vivid i, Vivid e, Vivid S5, Vivid S6, Logiq e, Logiq i, Logiq Book, and Logiq Book XP, the GE 3S-RS is a versatile ultrasound transducer. The GE 3S-RS can be used for cardiac applications including transcranial, abdominal, and cardiology. With a 1.7 – 4.0 MHz broadband multi-frequency range, and ergonomic design, the GE 3S-RS makes diagnostic imaging easier for users.

Features

- Wide-band phased array
- Adult and pediatric cardiac, abdomen, OB/ GYN, thoracic/pleural
- 19.3 x 27.6 mm footprint
- 1.7 - 4.0 MHz imaging frequency
- Multi angle

