



The **GE 4DE7C Transducer** is a wide band convex volume ultrasound transducer probe. The GE 4DE7C is used with OB/GYN and urology clinical applications and exams. This transducer has a frequency range of 3.5 to 10.0Mhz and is compatible with the Logiq series ultrasound systems. The GE 4DE7C transducer can be purchased from Soma Technology, Inc. fully refurbished.

## Features

- 3.5 – 10.0 MHz
- 32 x 27 mm Foot Print
- 116° Volume: 146° x 90° Field of View



# Specifications

## Physical Specifications

**Applications:** OB/GYN, Urology, Endocavitary

**Footprint:** 22 x 23 mm

**Bandwidth:** 4-9 MHz

**Biopsy Guide:** Single-angle, reusable

