



The **GE CARESCAPE ONE** is a patient monitor intended for use in multiple areas and intra-hospital transport within a professional healthcare facility. The CARESCAPE ONE multi-parameter monitor is indicated for the monitoring of hemodynamic (including ECG, ST segment, arrhythmia, detection, invasive pressure, non-invasive blood pressure, SpO₂, pulse rate, and temperature), and respiratory (impedance respiration and CO₂ airway gas) physiological parameters.

Features

- Adult, pediatric, and neonatal use.
- Can be used as a full-standalone monitor.
- Small and compact size support great intra-hospital transport.
- Touchscreen display.



