



The **Carestation 620** Anesthesia Machine is one of the newest additions to the GE anesthesia range. Its sleek compact design makes it ideal for operating rooms. The Carestation was designed specifically for low-flow anesthesia delivery and fast wash-in/wash-out requirements. Making it easy to ensure a quick patient response to changes in fresh gas or anesthetic agents.

The 620 comes with a 15 inch touchscreen display that features intuitive CARESCAPE inspired interface. Including a Bi-level work surface that offers illumination and intelligent flow control lights to clearly indicate active flow controls when in use.

Features

- Elegant modern design in a slim, compact frame well suited for constrained environments
- Simple and easy to use 15" touchscreen ventilator display
- Intuitive CARESCAPE™ inspired user interface for the unified Carestation user experience
- Integrated CARESCAPE Respiratory Module
- Time saving tools to help streamline clinician workload
- Scalable software and hardware features: "build your own" Carestation
- ecoFLOW display option may help clinicians mitigate the risk of hypoxic mixtures while helping to reduce agent use by using low and minimum flows with continuous gas monitoring



Specifications

Physical Specifications

Dimensions

Height: 135 cm / 53.1 in

Width: 82.5 cm / 32.4 in

Depth: 75 cm / 29.5 in

Weight: 145 kg / 320 lb

Top shelf

Weight limit: 25 kg / 55 lb

Width: 41.3 cm / 16.3 in

Depth: 27.0 cm / 10.6 in

Work surface

Height: 83.6 cm / 32.9 in

Size: 1930 cm² / 299 in²

Size: 2950 cm² / 471 in²

(with optional flip shelf)

Upper left Datex-Ohmeda (DO) dovetail

Dovetail length: 54 cm/21.3 in

Lower left Datex-Ohmeda (DO) dovetail

Dovetail length: 28 cm/11.0 in

Right Datex-Ohmeda (DO) dovetail

Dovetail length: 96.4 cm/38.0 in

Drawers (internal dimensions)

Height:

Top and middle: 8.6 cm/3.4 in

Bottom: 13.3 cm/5.2 in

Width: 34 cm/13 in

Depth: 37 cm/14.6 in

Manual ventilation bag arm (optional)

Arm length: 39.8 cm/15.7 in

Bag arm height (adjustable): 53 cm/20.9 in; 136 cm/53.5 in

Casters

Diameter: 12.5 cm/4.9 in

Brakes: Central Brake

