

GE Corometrics 259 Fetal Monitor



Features:

- GE-exclusive Smart BP* delays blood pressure measurement during contractions to reduce patient distraction and provide more meaningful diagnostic data.
- ChartLight with backlit keys and a dimmable display helps you work comfortably while your patients remain relaxed in a low-light environment.
- Watertight transducers have 9-crystal ultrasound technology that allows for easier pick-up and less repositioning, even for difficult to monitor patients.
- 9-crystal design achieves a wider focal region for more uniform coverage at greater depths than traditional 7-crystals.
- Minimizes monitor interactions



Typical Manufacturer's Image

GE Corometrics 259 Fetal Monitor

Specifications

| | |
|--|---|
| Maternal blood pressure: | DINAMAP* SuperSTAT* non-invasive blood pressure technology speeds up determinations by as much as 66% after the first reading to enhance patient comfort and allow you faster access to patient data. |
| Maternal temperature: | Integrated Exergen Temporal Arterial Thermometer allows you to obtain maternal temperature using non-invasive infrared technology. Readings automatically print on the fetal strip. |
| Maternal oxygen levels: | The monitor provides a choice of integrated Masimo® SET or Nellcor® maternal pulse oximetry, making it easy for you to monitor, document, and even view the waveform. |
| MECG: | Selectable 3 Lead maternal ECG allows clinicians to provides clinicians QRS wave amplitude and meaningful diagnostic data. |
| Additional technology includes: | <ul style="list-style-type: none">• Centricity* Perinatal (CPN) compatibility, which streamlines the capture and archiving of patient data and helps speed workflow and ensure data availability.• Telemetry Interface allows for continuous monitoring without changing transducers during ambulation or laboring in water. |



Soma Technology Inc.
166 Highland Park Drive
Bloomfield, CT 06002 • USA
1.800.438.7662
www.SomaTechnology.com

Note: The technical data given in this publication is for general information and are subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.