



This **Aladin 2 Cassette** delivers a isoflurane anesthetic agent electronically from the host anesthesia delivery system. The Aladin Cassette can be used on the GE Aisys CS2 anesthesia machine. Since these anesthesia vaporizers are electronically controlled, they feature automatic recordkeeping of set agent concentrations and agent usage calculations. The GE Aladin Cassette weighs in at less than 7.5 pounds and is easily carried. The isoflurane Aladin 2 cassette comes in 1 configuration related to refilling the anesthetic agent. It comes with the Easy-Fil.

## Features

- Prismatic site glass for clear indication of agent level.
- Automatic record-keeping on the anesthesia machine.
- Overfill protection.
- Electronic control.



## Specifications

### Dimensions

**Height:** 2.76 in. (7 cm)  
**Width:** 5.51 in. (14 cm)  
**Depth:** 9.45 in. (24 cm)  
**Weight** (empty): 6.2 lbs (2.8 kg)

### Agent Capacity

**Total:** 220 mL  
**Normal Fill** (when indicator shows empty): 125 mL

### Filling

**Overfilling Protection:** Overfilling prevention system built into the cassettes  
**Filling system:** Easy-Fil  
**Filling Speed:** > 2 mL/s

### Electronic Specifications

**Setting Ranges:** OFF, 0.2 to 5% in fresh gas flow, resolution 0.1%  
**Accuracy:** 18° to 25°C  $\pm$ 10% of setting or  $\pm$ 0.2% v/v (whichever is greater)  
Fresh gas flow range 1 to 10 L/min; 10° to 35°C  $\pm$ 20% of setting or  $\pm$ 0.4% v/v (whichever is greater) Fresh gas flow range 0.2 to 15 L/min

### Environmental

**Operating range:** 10° to 35°C/50° to 95°F  
**Storage range:** -25° to 50°C/-13° to 122°F

