



The **GE Datex Ohmeda Tec 4** anesthesia vaporizers are designed for “out of circuit” use in continuous flow techniques of inhalation anesthesia. The vaporizers are temperature and flow compensated so that their output remains relatively constant despite cooling due to vaporization and variations in inlet flow. A single control dial, with a concentration scale calibrated in % of anesthetic agent vapor per total volume (v/v).

Features

- Weight: 7.6 kg / 16.785 lbs
- Flow Range (LPM): 0.25-15
- Ambient Temperature Range: 18-35°C
- Capacity (dry wick): 135 cc
- Capacity (wet wick): 100

