



The **GE Logiq V2** is a laptop-style, portable color ultrasound system. The Logiq V2 was designed to fit all your ultrasound needs for a variety of applications. From Abdomen, Obstetrics, Gynecological, Vascular, Transcranial, Musculoskeletal, Urological, Small Parts, Pediatric, and Neonatal imagery. With color harmonic imaging the Logiq V2 can provide high-resolution images within deep penetration. Additionally offering real-time Speckle Reduction imaging that can reduce the noise and enhance the detail of the tissue for better diagnostic imagery. The ultrasound has a 15 in LCD high-resolution screen and an intuitive keyboard with programmable presets.

## Features

- Cart for additional stability with additional probe port and probe holder options
- Protective wheeled case for easy transport outside the office
- Wireless adapter for added portable connectivity
- On-board reporting to add key images to a report with exam findings and results that can be shared with an EMR
- MyTrainer on-board support program for real-time guidance and support



# Specifications

## **Easy-to-Use Color Imaging in a Portable Ultrasound System**

The compact LOGIQ V2 includes many advanced imaging technologies from decades of GE experience, including:

- Coded Harmonic Imaging enables high resolution even at deep penetration.
- Real Time Speckle Reduction Imaging (SRI-HD) reduces noise while enhancing true tissue detail.
- CrossXBeam™ enhances your ability to visualize tissue interfaces and differentiate borders.

## **Push Just a Few Buttons, Enhance Exam Efficiency**

With the LOGIQ V2, you get fast, accurate information with few keystrokes. For example:

- SonoBiometry – Performs key fetal biometry measurements in seconds.
- Auto IMT – Automates measurement of intima-media thickness for fast carotid exams.
- Scan Assistant – Puts imaging protocols at your fingertips to reduce keystrokes by as much as 22 percent
- Auto Optimization - Optimizes B-Mode image to help improve contrast resolution.

## **A System as Versatile as Your Practice**

The versatile LOGIQ V2 system accommodates a wide range of exams, including abdominal, obstetrics and gynecology, small parts, vascular, cardiac, pediatric and transcranial.

## **Scanning Support - Real Time, Real Easy**

The LOGIQ V2 gives you access to real-time guidance to support exam efficiency and accuracy:

- Scan Coach – Innovative scan plane reference tool that displays reference images, animations, and schematics during live scanning.
- Onboard help – Interactive access to information on settings, operation, connectivity, and maintenance, plus a link to the user manual for more details.

Plus, you become a member of the LOGIQ Club, with online access to tools and resources to help you experience the full power of your LOGIQ ultrasound system.

