



The **GE MAC 5** EKG machine is a portable electrocardiograph that can capture 3, 6, or 12-lead electrocardiograms on adult and pediatric patients. It features interpretive analysis and the ability to connect to a network using a wired LAN connection or wireless WiFi access point. The MAC 5 comes in 3 different models: the Mac 5 A4 with an A4 printer, the MAC 5 A5 with an A5 printer, and the MAC 5 Lite with no printer.

## Features

- 3, 6, and 12-lead ECG
- 8.9-inch touchscreen
- Gender and pediatric-specific interpretations
- Easy-to-clean design





