



The NuBOOM Ultra 2 or 4 is the most complete tower visualization product from CV Medical. With optional suspended lights and radiation shields, plus a touchscreen tablet control panel for the monitors, you can see and visualize multi-modality images with ease. Rapid renovation from procedure room to surgical theater with NuBOOM Ultra.

## Features

- High Definition LED displays
- Four 26" or two 24" LED monitors
- Simultaneously presented modality images on split screens – up to 8 images on 4 monitor unit
- Displays can be positioned at 7 points of articulation to provide comfortable 15 degree downward viewing angle
- Two equipment rails for mounting shelves and accessories
- Up to 15 image source inputs
- Two positionable booms support 4 display monitors
- Articulating boom arms surround up to a 360° surgical area
- Accommodates multiple monitor viewing from a variety of angles; Booms and arms retract for open surgeries
- Available with HD-SDI (High Definition Serial Digital Interface) input and HD display monitor configurations
- Optional wall displays
- Optional touch panel locations
- Upgradable and relocatable

