



GE OEC Apex CV pain management table/C-arm table with hand control

Features

- Efficient Low Attenuation Design
- Lighter Easy Transport
- Easy Maintenance
- Tilt, Roll and Float
- 30" Table Elevation
- 84" Head-to-Toe Imaging



Specifications

Dimensions

Table Top Length (Retracted): 90" (228 cm)
Table Top Width: 30" (81 cm)
Dimension of Image able Area: 18" x 84" (45 cm x 213 cm)
Weight: 850 lbs (385 kg)

Electrical

Power Requirements: 120V @60 Hz 10 Amps, 240V @ 50 Hz 5 Amps

Patient Load Capacity

400 lbs (180 kg)

Movements

Lateral Travel: +/-6" (15 cm)
Longitudinal Travel: 48" (122 cm)
Trendelenburg Travel: +/-20°
Lowest Table Height: 30" (76 cm)
Highest Table Height: 42" (106 cm)
Lateral Roll/Tilt: +/-20°

General

Table Top Composition: Carbon Fiber
Table Top Radiolucency: .55mm AL Equivalent
Table Pad Radiolucency: .50mm AL Equivalent (1/2 pad std.)
Approvals and Certifications: CE, VDE, CSA Approvals Pending

Accessories Included

1/2" Low Attenuation Table Pad
Single C.F. Arm Board with Pad
C.F. Tabletop Accessory Plate
Restraint Strap
Wristlets
Control Panel
Control Panel Covers (1 case of 25)
Hand Control

Accessories Available

Dual C.F. Arm Board with pad
1" Extra-Comfort Table Pad
Under-table mounted Radiation Shield I.V. Pole
Anesthesia Screen Holder

