



The GE Healthcare Venue Ultrasound is designed with a focus on emergency and critical care, as such health care providers have unique needs but often contend with ultrasound equipment from other medical disciplines. Specifically, Venue offers a Shock Toolkit of automated tools – developed on principles of machine learning – to gain information about the heart, lungs and inferior vena cava, automating the parts of shock evaluation that are the most tedious or difficult. The Auto VTI tool can help the user to visualize the patient’s trend and to determine the next course of action in treatment. The Venue features a seamless flat display; up to four-hour battery life; 19-inch articulating monitor; four probes with probe hangers and more.

Features

- Seamless display. For easy cleaning.
- You can focus on the patient. Simple and intuitive to use.
- ECG and TEE functional.
- Secured by a multi-year warranty. * Support you can count on.
- Suspension with large wheels and small footprint. Designed for robust use in emergency rooms. Just push the Venue where you need it.
- Handy key probe. Control the system with freely programmable buttons on the probe handle.
- Investigations that actually require three hands can be done with two hands.
- Ready to use when you need it. Up to 4 hours of battery power and a clearly visible remaining time display that changes color as the battery level gets lower.
- Large 19-inch LCD multi-touch monitor with height-adjustable swivel arm. You can position the monitor so that you always have the best view of it.
- Four probe connections. For the variety of examination possibilities.
- Intelligent cable management. Four connected and ready-to-use probes. Well-placed brackets help keep cables out of the ground.
- Network technology with state-of-the-art security standards. Whether wired or optional via WLAN - integrate the Venue easily into your network.



Specifications

Monitor

Instead of the stylus-based & rudimentary touchscreens of the 40 and 50, the Venue has a large, 19-inch articulating monitor that allows you to see images from across the room and reposition for various procedures. The flat-screen display is also seamlessly one piece of glass, making cleaning a breeze. There are also no buttons, which let you go right to work.

Cart

The new Venue cart is designed for cable management and efficiency. It has four probe ports in the back allowing you to have more transducers at your disposal. There are also integrated hangers that help keep the probe cables off the ground, helping them last longer. The cart also comes with cardiac-specific tools, like TEE and ECG. Its wheels have also been improved. They are bigger and easier to maneuver the Venue where you need it to go.