



The **GE Voluson S10** is an ultrasound system designed for basic imaging along with Obstetrics and gynecology. The Voluson S10 includes 2D, 3D, and 4D imaging with HD Live technology. This technology allows for a better rendering of images, creating an amazingly realistic image of a fetus through sonographic data. The System offers advanced tools that include Speckle Reduction imaging, CrossXBeamCRI, HD-Flow, and Dual-view. The GE Ultrasound provides archiving and reporting, making it easy to transfer patient data. This is done through a DICOM connection to an electronic medical record system or via a USB or DVD. Data and images can also be texted or emailed directly to the patient for convenience and patient satisfaction.

Features

- **More confidence:** Outstanding image quality, automated measurement tools and easy to use volume imaging.
- **More image information:** the large Voluson S10 monitor displays 98% more image information than previous Voluson Signature Series models.
- **More workflow speed:** the Voluson S10's intuitive system design helps reduce keystrokes and keeps critical key functions at your fingertips.
- **More scanning efficiency:** Four active probe ports provide easy access to ultra-light weight probes at the touch of a button.



Specifications

System Overview

Monitor (inch)	23" LED
Touch Screen (Inch)	No
Trackball or Trackpad	Trackball
CP Back-Lighting	Yes
Battery	Yes
Independent Steer & Lockable Wheels	Yes
Connecting	Text, email and USB

