



Increase your surgeons' operating room performance with the **Great White Maxima Stirrups** which permit simple adjustment of lithotomy and abduction intraoperatively. Exclusive axial rotation from the handle for superior patient positioning. Great White Maxima Stirrups range of motion is -35° low lithotomy to 84° high lithotomy. Adduction -9°, abduction 25° and also externally rotates out 20°, with the control to lock the stirrup anywhere in between. The Great White ships fully assembled. Boot pads are included. Two required SpringLoc Clamps sold separately.

## Features

- NEW FEATURE! Exclusive axial rotation from the handle for superior patient positioning
- A wide boot and a stronger gas cylinder make this stirrup perfect for your heaviest patients
- 800 lb (363 kg) surgical patient weight capacity
- Boot free-floats throughout each range reducing the risk of superficial nerve injury
- Lithotomy and length indicators help to position and reposition the legs
- Extended lateral boot design helps reduce the risk of peroneal nerve injury

