



High quality, feature-rich, and priced right, the **Great White Premium Stirrups** are exactly what your surgical team requires. The enhanced range of motion allows for -35° for low lithotomy, to 84° high lithotomy; abduction -9° to abduction 25° and also externally rotates out 20° . Two required SpringLoc Clamps sold separately.

Features

- NEW FEATURE! Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -35° low to 85° high; Adduction -9° and abduction $+25^{\circ}$
- 400 lb. surgical patient weight capacity
- Included Premier boot liners wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- Lithotomy indicators help to accurately reposition legs

