



The **HP Codemaster XL+** is a high performance Defibrillator monitor. The HP Codemaster XL+ offers multi-functional electrodes that allow you to do pacing and defibrillation with the same pads. At the same time, the system offers a hands-free operation that allows for remote defibrillation. The Defibrillation system has built-in paddles for pediatric and adult patients. With its easy 1-2-3 operation, the system is easy to use and offers a speed response in an emergency. The system can charge from 0 to 360 in 5 seconds and store up to 28 events in the memory. The HP Codemaster XL Plus has a bright 5-inch screen that can tell you the condition of the patient at a glance.

Features

- Defibrillator with ECG and recorder
- Provides 1-2-3 operation
- Built-in pediatric paddles under slide-off adult paddles
- Preset heart rate alarms
- High-capacity battery
- Built-in diagnostics



