



The Haemonetics Cell Saver Elite+ is an autotransfusion device that allows medical facilities to collect and salvage the patients blood during and after a procedure to then deliver it back to the patient after it has been processed and washed. Autotransfusion systems are typically used in procedures where high blood loss occurs. Soma Tech Intl offers refurbished Haemonetics Cell Saver Elite+ devices. Soma refurbishes cell savers in their own ISO:13485 certified facility and is not affiliated with Haemonetics.

Features

- Reduces the need for allogeneic transfusions with autologous transfusion.
- Packs and washes red blood cells with over a 90% recovery rate.
- Touch-screen color monitor.
- Customizable settings.



Specifications

Dimensions

Height: 72" (183 cm)
Width: 21" (53.3 cm)
Depth: 26.5" (67.3 cm)
Weight: 95 lbs. (43.1 kg)

Power

Power Supply: 200 - 240 VAC
Fuse: T 3.15 A / 250 V
Operating frequency: 050 - 60 Hz

