



The Hemochron 801 is a portable, dual-well whole blood coagulation monitoring instrument. The 801 provides on-screen plasma equivalent PT and APTT test results. Portable, battery-operated coagulation instruments designed for bedside use.

Features

- Hemochron 801: Dual well analyzer
- Portable
- Battery-operated coagulation instruments designed for bedside use
- 2 Channels
- Range of Coagulation Test: APTT, PT, TT, HNTT, ACT, PDA



Specifications

Dimensions

Height: 15 cm
Width: 25 cm
Depth: 23 cm
Weight: 2.66 kg

Number of Channels

1 Channel

Timing Range

0 - 1500 Seconds

Incubation Temperature

37°C ± 0.5 °C

Incubation Warm-Up Time

4 - 6 hours

Operating Time on Full Battery Charge

4 - 6 hours

Line Voltage

120 volts AC % Hz)
100, 220, 230 optional

Power

11 watts

Range of Coagulation Tests

APTT, PT, TT, HNTT, ACT, PDA

