



The HEMOCHRON® Model 401 is a portable, single-well whole blood coagulation monitoring instrument. The 401 provides on-screen plasma equivalent PT and APTT test results.

## Features

- Fully portable system weighing only 1.7 kg
- Uses fresh whole blood
- Test results within minutes
- Timing range is 0 to 1500 seconds
- Operates on AC current or self-contained rechargeable batteries

## Test Menu:

- Celite® Activated Clotting Time (ACT)
- Kaolin ACT
- Glass Activated ACT
- Activated Partial Thromboplastin Time (APTT)
- Prothrombin Time (PT)
- Thrombin Time (TT)
- Heparin Neutralized Thrombin Time (HNNTT)
- High Dose Thrombin Time (HiTT)
- Fibrinogen
- PDAO (Protamine Dose Assay)

