



### Models Available

- Hemochron 801: Dual well analyzer
- Hemochron 401: Single well analyzer

Portable, battery-operated coagulation instruments designed for bedside use.

### Specifications:

Hemochron Model	401	801
Number of Channels	1	2
Timing Range	0 – 1500	0 - 1500
Incubation Temperature (°C)	37 ± 0.5	37 ± 0.5
Operating Time (hours)	4 - 6	4 - 6
Line Voltage	120	120
Dimensions (w x d x h, cm)	18 x 23 x 12.5	25 x 28 x 15
Weight (kg)	1.7	2.6
Range of Coagulation Tests	APTT, PT, TT, HNTT, ACT, PDA	APTT, PT, TT, HNTT, ACT, PDA

