



The **Hill-Rom VersaCare** is a versatile med-surg bed. The VersaCare bed is a low-high bed that can fully be adjusted to fit the patient's needs. The Med-Surg bed includes a one-button control system. This makes it easy for the patient to adjust the bed into any position for maximum comfort.

Features

- One -Button controls help the patient move the bed into positions
- 3 Level Bed alarm - can sense if the patient is lying in bed, sitting on the edge, or is out of the bed
- Adjustable length of the bed to help prevent patients from sliding down the bed
- Low Chair Function allowing the patient to be placed in a sitting position while the bed is in a low position
- Smart Bed technology that can connect to the nurse call system or EMRs at your facility
- Headboard and push handles are stationary in a position while the patient can be adjusted
- Hands-free foot controls for caregivers to adjust the patient
- Integrated battery backup for when patients are being transferred
- CPR release lever can instantly lower the head of the bed for chest compressions
- The bed can easily be moved into 5 different functional positions
 - High/Low
 - Fowler
 - Knee
 - Trendelenburg
 - Reverse Trendelenburg



Specifications

Bed Height

Low position – floor to the center of the deck: 13 in (33 cm)
 Low position – floor to the edge of the deck: 17 in (43.2 cm)
 High position – floor to the center of the deck: 30 in (76.2 cm)
 High position – floor to the edge of the deck: 34 in (86.4 cm)

Sleep Surface

Width: 35.5 in (90.2 cm)
 Length, fully extended: 86 in (218.0 cm)
 Length, fully retracted: 75 in (191.0 cm)
 Radiolucent sleep deck: No
 Surface weight limit: 500 lb

Overall Width

Siderails stowed: 37 in (94.0 cm)
 Siderails up: 40 in (102.0 cm)

Overall Length

Fully extended: 94.5 in (240.0 cm)
 Fully retracted: 82.5 in (210.0 cm)

Bed Angle

Head section: 0° to 65°
 Knee section: 0° to 16°
 Foot section: 0° to -27°
 Max Trend/Reverse Trend: 15° / 10°
 Head and Trend angle indicators: Yes

Casters and Breaking System

Optional IntelliDrive power transport: Yes
 Caster diameter: 5 in (12.7 cm)
 with IntelliDrive power transport: 6 in (15.2 cm)
 Four-wheel braking: Yes
 Brakes lock both rolling and swiveling: Yes

Patient Weight Limit

Patient weight, all surfaces: 500 lb (227 kg)
 Safe working load, including all accessories: 550 lb (250 kg)

