



Features

- Low height improves on/off transfer and increases patient safety.
- One-Step® TuckAway siderails for zero-gap transfers.
- Ergonomic push handles to reduce caregiver strain during transport.
- GentleRide™ shock absorbing frame and siderail system for improved patient comfort



Specifications

Overall length – standard	83" (210.8 cm)
Overall length – surgical	92" (233.7 cm)
Overall length – OB/GYN	80.5" (204.5 cm)
Overall width (siderails up) – standard	32" (81.3 cm)
Overall width (siderails up) – wide/electric	36" (91.4 cm)
Overall width (siderails up) – OB/GYN	35.5" (90.2 cm)
Upper frame width (siderails stored) – standard	26.5" (67.3 cm)
Upper frame width (siderails stored) – wide/electric	30.25" (76.8 cm)
Upper frame width (siderails stored) – OB/GYN	31" (78.7 cm)
Overall siderail length – standard	47" (119.4 cm)
Overall siderail length – OB/GYN	37" (94 cm)
Siderail height above sleep surface – standard	14" (35.6 cm)
Siderail height above sleep surface – OB/GYN	13" (33 cm)
Mattress size – standard	26" x 75" (76.2 cm x 190.5cm)
Mattress size – wide/electric	30" x 75" (76.2 cm x 190.5cm)
Mattress size - surgical	26"x78"(66cmx198.1cm)
Mattress size -OB/GYN	30" x 76.5" (76.2 cm x 194.3cm)
Mattress thickness	3" (7.6 cm), 4" (10.2 cm), or 5" (12.7 cm)
Hydraulic cylinder specifications: Rated pressure	1000 psi (6895 kPa)
Hydraulic cylinder specifications: Burst pressure	5000 psi (34474 kPa)