



The Affinity Four birthing bed delivers on the promise of a safe birthing bed that functions in the busiest environment with ease. With the Affinity Three Stow and Go™ feature, the integrated foot section can simply slide under the bed during birth and back into place afterwards. Or, if you choose the ergonomically designed lift-off version, simply place the foot section on the built-in stand.

Features

- The Stow and Go foot section slides and stores under the bed – no lifting required. The foot section easily slides back into place without false latching or interference with linens.
- With battery backup, caregivers can use all the bed's features, including changing height and position, when away from a power outlet, reducing the time from "decision to incision" in critical situations.
- The gas-assisted Trendelenburg mechanism allows head of bed to immediately lower to 8 degrees.
- OneStep siderail release provides zero transfer gap for safer and easier patient transfers.
- The bed's streamlined design, with all the equipment stored completely within the perimeter (instead of protruding out), protects caregivers as they move quickly around the patient.



Specifications

Physical Description

Total Length: 90" (229 cm)
Length from roller bumpers to break in seat section: 61 7/8" (157 cm)
Length from break in the seat section to end of bed: 28 1/8" (71 cm)
Maximum Width (siderails stored): 36" (91 cm)
Maximum Width (siderails up): 39" (99 cm)
Maximum Headboard Height: 52" (132 cm)
Maximum Siderail Height (without mattress): 14 3/4" (37.47 cm)
Minimum Under-Bed Clearance: 5" (13 cm), 3" 7.6 cm for Stow and Go™ Beds
Wheel Base: 50" x 29" (127 cm x 74 cm)
Mattress Width: 35" (89 cm)
Mattress Length: 78" (198 cm)
Mattress Thickness (head/seat): 4" (10 cm)
Mattress Thickness (foot): 3" (8 cm)
Detachable Power Cord, IEC 320/interface: US, 84" (213 cm)/international, 98" (249 cm)
Caster Size: 6" (15 cm) or 8" (20 cm)
Total Weight: 480 lb (218 kg)
Head Section Inclination (maximum): 65°
Seat Section Inclination (maximum): 15°
Bed Height Range 18" (46 cm) to 34" (86 cm)
Bed Height Range (with mattress) 22" (56 cm) to 38" (97 cm)
Trend-Like position (maximum): 8°

Environmental

Use:
Temperature 59°F (15°C) to 104°F (40°C) ambient temperature
Relative humidity range 10% to 85% non-condensing
Transport & Storage:
Temperature -40°F (-40°C) to 158°F (70°C)
Relative humidity 95% non-condensing
Pressure 15 "Hg (50 kPa) to 31 "Hg (106 kPa)

Fuse

Air system fuse (air system optional): 2 A, 250 V~, 5 x 20 mm, UL 198G Fast Acting
Battery fuse 10 A, 32 V~, ATO
Mains fuses (120 V bed model) 4 A, 125 V~, 5 x 20 mm, UL 198G Slo-Blo@a or equivalent
Mains fuses (230 V bed model) 2 A, 250 V~, 5 x 20 mm, IEC127 Sheet III, Time Delay