



The **Hillrom Advance Series Beds** are 4-motor retractable type beds that move the patient closer to the head end of the bed automatically when the head section is elevated. The bed is equipped with a night light that automatically illuminates an area under the bed, at both sides, when the ambient light decreases. The Advance Series beds are mechanically built with an automatic contour position to make the patient feel more comfortable.

Features

- 4 individual motors to power separate functions of the bed.
- Maximum patient weight of 500 lb.
- Automatic Contour.
- Emergency manual operation.
- 4 precision bearing, swivel-type casters.
- Hilow, Automatic Contour, Chair Position, and Trend./Reverse Trend. movements.



Specifications

Dimensions

Max Height: 43" (109.2 cm)

Width: 41.5" (1105.4 cm)

Length: 93" (237.5 cm)

Max Patient Weight: 500 lbs (227 kg)

Movements

Head Section Inclination

Knee Section Inclination

Automatic Contour Up

Automatic Contour Down

Automatic Contour

Hilow Range

Trendelenburg

Reverse Trendelenburg

Nurse Control Panel

Trendelenburg/Reverse Trendelenburg

Chair Positions

Hilow Control Switches

Prevention Selector Switch (Optional)

Comfort Selector Switch (Optional)

Mattress Pressure Adjustment Switches

Pressure Relief Selector Switch

Heel Relief Selector Switch

Siderail Back Lighting

Siderail Hilow Function

Siderail Knee Function

Scale

Siderail Sleep Surface Lockout

Sleep Surface Power Lockout

Bed Motor Power Lockout

