



The **Hillrom Advanta P1600A/B** are four-motor, retractable-type beds with multiple articulation features. A variety of options are available, including nurse call, patient scale, a patient position monitoring system, night light, air system, and side rail controls. This retractable-type bed means that the articulating frame moves towards the head end of the bed as the head section is elevated. This has the effect of keeping the patient in the same relative position with regard to objects outside the bed, such as shelves or night stands.

Features

- 4 individual motors to power separate functions of the bed.
- Articulating Frame.
- Automatic Contour.
- Emergency manual operation.
- 4 precision bearing, swivel-type casters.
- Hilow, Automatic Contour, Chair Position, and Trend./Reverse Trend. movements.



Specifications

Dimensions

Max Height: 48" (122 cm)
Width: 42.75" (108.59 cm)
Length: 91" (231 cm)
Max Weight: 500 lbs (227 kg)

Movements

Head Section Inclination (maximum): 65°
Knee Section Inclination (maximum): 35°
Automatic Contour Up (maximum): 15°
Automatic Contour Down: 28°
Automatic Contour: 10°
Hilow Range: 13.25" (33.66 cm)
Trendelenburg: 11°
Reverse Trendelenburg: 11°

Patient Siderail Controls

Standard

Head up/Head down (with automatic contour)
Knee up/Knee down

Optional

Nurse call
Nurse call/Answer indicators
Room/Read lights
TV on/off
Channel up/down
Radio on/off
Select (closed captioned TV)
Volume up/down
Speakers

Caregiver Siderail Controls

Standard

Bed up/Bed down
Head up/Head down (no automatic contour)
Knee up/Knee down
Lockout indicator
Head angle indicator

Optional

Nurse call
PPM mode: position, exiting, out of bed
PPM volume: low, med, high
PPM enable control
PPMmode indicators: position, exiting, out of bed
Brake not set indicator
Bed not down indicator

