



The **Hillrom P8005** transport stretcher is used to transport patients in a variety of different medical facilities such as hospitals and surgical centers. The Hillrom P8005 stretcher is easy to steer and stop for maximum control. It features a low height for safe ingress and egress, and pedals for hilow, trendelenburg, and reverse trendelenburg movements. This Hillrom transport stretcher has a maximum patient weight capacity of 700-pounds.

Features

- Dual locking 4-wheel brakes.
- Siderails that fold down.
- Head and foot brake and steer pedals.
- Standard 3-inch mattress.



Specifications

Dimensions

Height: 23" (58.4 cm) - 34.25" (87 cm)

Width: 32" (81.3 cm)

Length: 84.5" (214.6 cm)

Weight: 265 lbs (120.2 kg)

Standard Mattress

Height: 3" (7.6 cm)

Width: 26" (66 cm)

Length: 75" (190.5 cm)

Weight: 13 lbs (5.9 kg)

Movements

Height: 23" (58.4 cm) - 34.25" (87 cm)

Max. Head Elevation: 90°

Max. Trend/Reverse Trend: 18°

Min. Trend/Reverse Trend: 12°

