



The **idi Aspect 100U** is a multi-purpose mobile urology table with motorized elevation, trendelenburg, tilt, and longitudinal motions. The Aspect 100U provides a large radiolucent imaging area for easy C-arm access and greater anatomical coverage. Urological accessories can easily be removed, and a lightweight radiolucent tabletop extension can be equipped to convert the table to multi-purpose imaging or general patient positioning.

Features

- Low X-ray attenuation radiolucent tabletop
- 2-inch thick memory foam pad
- Foot switch and hand pendant controls
- Standard US side-rails for different accessories and extensions



Specifications

Dimensions

Height: 31.3 - 41.1" (79.4 - 104.4 cm)
Width: 28" (71 cm)
Length: 46" (117 cm)
Weight: 686 lb (311.25 kg) (With extension)

Tabletop Extension

Type: Radiolucent
Width: 24" (61 cm)
Length: 34" (86 cm)

Movements

Controls: Foot switch and hand controls
Height: 31.3 - 41.1 in (79.4 - 104.4 cm)
Trendelenburg Tilt: $\pm 15^\circ$
Longitudinal Tabletop Travel: 10" (25.4 cm)

Imaging Area

Main Tabletop: 25"W x 20"L (63.5 x 50.8 cm)
Tabletop Extension: 24"W x 34"L (61 x 86 cm)
Total Imaging Area With Extension: 54" (137 cm)

Patient Weight Capacity

With Patient Upper Body on Main Tabletop: 525 lbs (238 kg)
With Patient Upper Body on Extension: 400 lb (181 kg)
Load Capacity on Extension Only: 200 lb (90.7 kg)

Electrical

100 VAC
60 Hz Operation (230 VAC optional)
 ≤ 10 amp @110V

