



The **Impact Eagle Uni-Vent 750** is a portable ventilator. It uses an electronically controlled microprocessor, allowing the patient's airways to be continuously monitored. The Uni-Vent was designed to meet the needs of transport and critical care. The System features a comprehensive alarm system and AWR (Air Worthiness Release), a standard for winged aircraft. The ventilator has multiple modes for ventilation including ACV, SIMV, and CPAP, and can be used on patients ranging from infants to adults.

Features

- Protective Soft Field Case
- AC Power Supply
- DC Power Supply
- O2 Hose, Green
- Operator's Guide

