



Soma offers the **Integra CUSA Excel+** refurbished up to 50% below the cost of new. The Integra CUSA Excel+ is an ultrasonic surgical aspirator that can perform fragmentation, irrigation, and aspiration simultaneously. It provides 23,00 to 36,000 cycles per second to the attached handpiece. The CUSA Excel+ features a closed recirculating water cooling system to reduce the heat generated by the high-frequency handpiece. The hollow surgical tip pushes sterile irrigation fluid through pre-aspiration holes and offers aspiration at the tip to remove irrigation fluid and tissue.

Features

- Fragmentation, irrigation, and aspiration.
- Automatically adjusts power to maintain proper tip vibration.
- User-friendly control panel.



Specifications

Dimensions

Height: 52 in (132 cm)
Width: 22 in (55.9 cm)
Depth: 28 in (71.1 cm)
Weight: 200 lbs. (93 kg)

Power

Voltage: 220-240V
Power: 110-120 VA

Liquid

Water Unmineralized water for cooling system
Vacuum On board pump

Cooling System

The device is equipped with a closed water cooling system

Environmental

Functioning temperature 10 to 40°C
Hygrometry 15 to 80% of humidity, without condensation

