



Soma Tech Intl offers the **Integra DUO** surgical headlight refurbished up to 50% below the cost of new. The Integra DUO is an LED headlight that provides optimal head comfort and brightness consistency. It can be used with a portable battery pack or be connected to AC power. Featuring an active cooling system, the onboard fan draws air through the channels through a duct system and exhausted behind the surgeon to cool off the LED lamp.

Features

- Ergonomic headlight with dual-layer padding and multiple adjustments for maximum comfort.
- A fully charged battery will operate for about 4 hours.
- The color temperature is 5900k.
- LED light source.



Specifications

General

Classification: FDA Class II, Exempt: CA, EU, AU: Class 1

Type: LED

Wattage: 9.3 Watts

Lamp Life: 48,000 Hours

Color Temperature: 5900K

Spot Size Range: 4.2cm to 12.5cm at 16 in (~40cm) working distance

Power Input: 9-15 VDC

Weight: ~14.9 oz (422.4 grams)

Battery

Type: Lithium Ion, 4-cell pack

Capacity: 48.2 Wh

Voltage: 7.2 VDC @ 6700mAh

Electrical Safety: IEC 62133 Ed. 2

