



Soma offers the **Integra MAYFIELD A2009 Ultra 360 Patient Positioning System** refurbished up to 50% below the cost of new. The MAYFIELD A2009 Patient Positioning System is where both the headrest and the skull clamp are attached to create a headrest system on a surgical table. It offers extra stability in keeping the patient in a fixed position for procedures that include Cranial, Spinal, and Neurosurgical. The positioning system has two handles that can rotate and self-locking independently for more flexibility for making adjustments. The patient positioning system A2009 is compatible with other MAYFIELD Skeletal stabilization Headrests, and Skull Clamps to create the unit you need for a procedure.

Features

- Handles that can rotate and are self-locking
- More movement and flexibility for quicker positioning with the skull clamps
- Reduced set-up time

