



Soma offers the Integra MAYFIELD XR2 A2079 Radiolucent Base Unit refurbished up to 50% below the cost of new. The MAYFIELD A2079 XR2 Radiolucent Base Unit can connect MAYFIELD Skull Clamps and Headrest to an operating table. Which is used for procedures like Cranial, Neurosurgical, or even spinal surgeries. The A2079 base is Radiolucent, making it compatible to be used in CT imaging and MR imaging needs. The Base Unit is compatible with other MAYFIELD Infinity Skeletal stabilization pieces.

Features

- Composed of carbon fiber to not interfere in the fluorography imaging
- East to set up and position with linking components
- Color-coded knobs on the unit to make adjustments and save time in the operating room
- No tools are needed to make adjustments
- Horizontal bar provides a wide imaging window-full, side-to-side lateral positioning (21.5-inches/54.6cm)

