



Soma offers the **Integra MAYFIELD A2010 Adult Radiolucent Horseshoe Headrest** refurbished up to 50% below the cost of new. The MAYFIELD A2010 Radiolucent Headrest is an adult headrest that can be used for Cranial, Spinal, and Neurological procedures. Providing non-invasive support while holding the head and neck of the patient in place in either the prone or supine positions. The headrest is made out of Carbon Fiber Radiolucent material so it can be used in CT and MR Imaging while producing minimal artifacts in the image.

Features

- Support for either prone or supine procedures
- Comfortably support that can be adjusted to fit the patient's head
- Vertical and Lateral Adjustments can be made to better position the patient
- Optional gel pads that are not included in the standard package.

