



Soma offers the **Integra MAYFIELD Infinity A1114 Skull Clamp** refurbished up to 50% below the cost of new. The MAYFIELD Infinity Skull Clamp holds the head and neck of the patient in a fixed position for either Cranial, Spinal or Neurosurgical procedures. The Infinity style of head clamp works in conjunction with the Infinity support system, along with a Horseshoe headrest. The head clamp sits slightly lower than the position of the headrest for better support of the patient's head during procedures.

Features

- Works with the MAYFIELD Infinity Support System
- Offering a lower profile, making it sit lower than the headrest
- Comparable with any of the Infinity Support System Horseshoe headrests
- Swivel arm that can rotate 360 degrees under pressure
- Offers multiple rocker arm positions: Standard, Triad, Pediatric, and Single-pin

