



Soma offers the **Integra MAYFIELD Infinity XR2 A2114 Radiolucent Skull Clamp** refurbished up to 50% below the cost of new. The MAYFIELD InfinityXR2 A2114 Radiolucent Skull Clamp holds the head and neck of the patient in a fixed position for either Cranial, Spinal or Neurosurgical procedures. The Skull clamp comes with 3 points for fixation of the patient's head, allowing for minimal artifacts. The Radiolucent Skull clamp comes with a focusing knob on the side. The XR2 A2114 works with a variety of MAYFIELD Infinity skeletal stabilization pieces to create the unit you need for a procedure.

Features

- Provides 3-point skeletal fixation to hold the head and neck in place
- Comparable to be used in Fluoroscopy, CT Imaging, and MR imaging areas
- MR imaging comparable up to 3 Tesla
- The extension is fixed by a double pawl, better control and less pressure needed to apply to the patients head

