



Soma offers the **Integra MAYFIELD Modified A1059 Skull Clamp** refurbished up to 50% below the cost of new. The MAYFIELD A1059 Skull Clamp offer 3-point skeletal fixation, holding the head and neck in a specific position for Cranial, Spinal or Neurosurgical procedures. Offering Flexibility in the positioning of the skull pins. The skull clamp sits slightly lower than the position of the headrest for better support of the patient.

Features

- Flexibility of the positioning of the skull pins
- 3-point skeletal fixation
- Offer a 80 lbs torque screw
- Quick release lock

