



Soma offers the **Integra MAYFIELD A2600M Spine Table Adaptor** refurbished up to 50% below the cost of new. The Mayfield A2600M Spine Table Adaptor can connect any of the Mayfield swivel adaptors, Skull clamps, or horseshoe headrests to a spinal table in preparation for a cranial procedure. Spinal tables like the Mizuho OSI Jackson Table are compatible and easy to set up through locking knobs that can connect through the rails brackets of the table.

Features

- Easy to connect to a spinal table
- Uses the locking knobs of the adaptor to connect to the rail brackets of the spine table
- Connects Mayfield and other compatible cranial stabilization tools
- Comparable pieces include Swivel Adaptors, Skull Clamps, and Horseshoe headrests

