



Soma offers the **Integra MAYFIELD A2101 Ultra Base Unit** refurbished up to 50% below the cost of new. The MAYFIELD Base Unit is a device that can connect to MAYFIELD Skull clamps and headrests to create the ideal combination for any Cranial, Neurosurgical, or spinal surgery. For patient positioning in prone, supine, or even a sitting position the Ultra Base Unit allows for heavier weight and more face-bearing. This improves the safety of the patient during procedures. The base unit offers width and locking positions that can be adjusted by using a wrench.

Features

- Increased weight and force bearing for patient safety.
- The Base tubing (15.25 inches) can provide more flexibility in lateral patient positioning
- Adjustable width of the end brackets from 6.75 inches to 8.25 inches can fit most tables
- Designed to be used with the MAYFIELD Tri-Star Swivel Adaptor (A2008)
- Can be used for procedures that require a prone, supine, or a sitting position

