



The **Leica M680 OH-2** Surgical Microscope has two independent zooms and two internal focusing systems, making the M680 ideal for co-operations. The M680 has a 2x2 formula that brings more individual freedom within the teamwork. Each surgeon can use a hand or foot switch to focus the image and well as engage either a high magnification for detailed observation or a low magnification for an overall view. The Leica M680 offers opportunities in reconstructive, hand, heart and spinal surgery, along with urology and gynecology.

Features

- Free floating Electronic Clutch-Release Floorstand
- M680 Optics
- Variable Optics 207mm-407mm
- Face to Face Bridge
- Two Binocular Inclinable 0-180
- Motorized Focus and Zoom
- 4 10X Widefield

