



The **Leica OH-1** surgical microscope is used for Neurological and spinal surgeries. Leica Microscopes are widely recognized for their optical precision and innovative technology. The OH-1 microscope has a Neurosurgical OP-Microscope that has fully motorized adjustment of the zoom and focus. The Binocular, swiveling, tube working distance 272-312mm, field of sight 20-72mm, cold light source 12V with 50W and $f = 325$.

Features

- OH1 Floorstand
- Inclinable 0-180° Binoculars
- 10X/22B Eyepieces
- Motorized Focus and Zoom
- New Fiberoptic Illumination
- Multifoc System
- Pistol Control

