



Features

- Used for angiography, CT /MRI /DSA scan and enhanced angiography in our body.
- High quality medical grade Polypropylene material
- Kit for CT injector systems
- Perfectly fit to CT injectors such as Medrad 200/300 vistron
- Standard pressure and highpressu for choosing.

General Specifications

- Applicable Scope: Liebel-flarsheim(USA)
- Optistar LE,Elite MRI High Pressure Injector
- Standard Configuration:2-60ml syringe,1-Y(T) connecting tube,2-suction tube (device)

