



Linvatec Universal Drive Console is a microprocessor-controlled power unit that drives and controls Linvatec handpieces. Console functions include: audio response feedback, normal on/off and variable (accelerator) footswitch speed controls, footswitch sensitivity control, handpiece aspirate mode, handpiece reverse/forward/oscillate cutting modes, RPM speed controls, and digital front panel display. A help menu provides functional explanations and surgeon custom programming. This drive console will accommodate large joint handpieces. 320W max. power. Mains 90-240V, 50/60Hz, 6A 456W

