



The **Maquet Alphaclassic Pro** mobile operating table is a general surgical table that can be used for cardiovascular, critical care or even surgical workplaces. The Alphaclassic Pro can lower to 23.6 inches and has a maximum weight capacity of 551 lbs. It has a longitudinal shift of 12.2 inches and a lateral tilt of 15° to either side.

Features

- Guide rails for inserting X-ray cassettes.
- 551 lbs. table weight capacity.
- Backlit hand control.
- Up to 5 subdivisions.



Specifications

Dimensions

Height: 23.6-37.4 in (60-95 cm)

Width: 23.2 in (59 cm)

Length: 84.3 in (214 cm)

Weight Capacity: 551 lbs. (250 kg)

Movements

Height: 23.6-37.4 in (60-95 cm)

Longitudinal Shift: 12.2 in (31 cm)

Table Width (across side rails): 23.2 in (59 cm)

Trendelenburg / Reverse Trendelenburg: 25°

Lateral Tilt to Either Side: 15°

Lower Back Plate: +70° / -40°

