



With a few positioning adjustments and optional accessories, the Maquet Alphamaxx 1133 is ideal for operations in gynecology, urology, orthopedics, traumatology and even neurosurgery. The table features electro-hydraulic actuators that control longitudinal shift, lateral tilt, and various table angles. The Alphamaxx 1133 adjusts from 23.4 to 41.6 inches in height, allowing clinicians to work comfortably from a sitting or standing position.

An easy click mounting system allows for a wide range of accessories to be quickly added to the table, efficiently transforming the chair for each unique procedure. The durable table's maximum load capacity is rated at 990 pounds.

Features

- 6-section table features head plate, upper back plate, lower back plate, seat plate and leg plates Unique split leg design for a wide range of positioning capabilities.
- Electro-hydraulically driven actuators control table movements.
- The motorized longitudinal shift for C-arm access.
- Individual leg plate adjustment, back, and headrest manipulation.
- Height adjustment allows for patient and clinician ergonomic comfort.
- Easy click mounting for a wide range of table accessories.
- The powerful table has a load capacity of 990 pounds.
- Support arms, leg clamps, joint cover, column casing and side rails made of durable, chrome-nickel steel.
- Corded hand remote, wireless controls, and footswitch available for the Maquet 1133.
- Emergency override panel on table column.
- Rechargeable battery and mains operation.
- 4 swivel casters for easy maneuvering.
- Mattress features 3.1-inch, SFC padding.



Specifications

Table Height	Raised: 41 1/2 Inches Lowered: 23 3/8 Inches
Trendelenburg	30° Reverse Trendelenburg: 30°
Head Section	Raise: 0-45° Lower: 0-45°
Lateral Tilt	0+-20°
Back Section	Down: 40° Up: 80°
Seat Section	Flex: 40° Reflex: 85°
Leg Section (Removable)	Down: 90° Up: 10°
Table Top Slide	9 inches
Return Level	'0' Button
Auto Drive	Hand Control
Brekes (Electric)	Hand Control Up: 10°
Weight Limit	990 lbs
Width	20 15/16 Inches With Rail: 22 3/4 Inches

