



The Maquet Lyra offers several patient positioning possibilities along with a wide range of angles and lateral tilt combinations. The Lyra Mobile OR Table was designed to work with almost all surgical disciplines and is compatible with existing Maquet Mobil Or table accessories that your facility may already have. The robust design of the table ensures the stability and safety for patients. The Maquet era has an overall load capacity of 360 kg (794 lbs) at the central position and 180 kg (397 lbs) without any restrictions.

Features

- Maximum Overall Load - It has a overall load (patient plus accessory) of 360 kg (794 lbs) in normal position, and 180 kg (397 lbs) without any restrictions
- Versatile angle, Tilt and heights - A 350 mm height range supports a variety of procedures, whether the surgeon is seated or standing. This improves ergonomics reducing strain and minimizing the risk of workplace injury.
- Design for universal use - A flexible mobile table can serve many purposes. Manual longitudinal shift improves intra-operative diagnostics and treatment. It delivers increased access for the surgeon while simultaneously improving accessibility for the C-arm. 360° imaging is possible when paired with carbon fiber accessories.
- Improved throughput between procedures - Gas-spring assisted leg plates can fold around the column, remaining attached during gynecological procedures. This reduces reconfiguration time and helps to increase throughput.
- High maneuverable - The special design of the castors allows the table to be used immediately, without being lowered onto a baseplate. The free-wheeling castors give the flexibility to move in both straight lines or laterally.
- Excellent clinical access - A concave base and smart design of the column and table top allow the surgical team to get close to the patient – without positional strain. The team can work comfortably and confidently.
- Safety in an Emergency - Maquet Lyra can still be used in case of a power failure or table software blackout. A hard-wired override panel is connected to a separate circuit for safe control. Table brakes can be manually released if necessary.



Specifications

Technical Data

Length of operating table top with standard back plate, head rest, and leg plate 2030 mm/80"

Length of operating table top without accessories 745 mm/29"

Width without side rails 540 mm/21"

Width across side rails 590 mm/23"

Net weight (without accessories) 190 kg/419 lbs

Maximum overall load (patient and accessories)

360 kg/794 lbs with restrictions

180 kg/397 lbs without restrictions

Adjustment Range

Height without padding 600–950 mm/ 24–37"

Trendelenburg/Reverse Trendelenburg 25°/25°

Lateral tilt 15° Lower back plate +75°/-40°

Leg plate (manual) -90° (down)

Manual longitudinal shift 310 mm/12"