



The **Medicap USB 200** is the medical industry's easiest full-featured video recording solution. Your videos and still images are recorded directly to a USB flash drive or optional hard drive. The recorded videos are DVD-quality, but without all the hassles of DVD disks.

Medical Video Recorder

Records from your endoscope, arthroscope, Surgical Camera, Ultrasound, etc. Saves to a flash drive or hard drive that can be accessed from any computer on your network.

Intuitive Use

It's so easy, no training is needed. Just plug it in and go. You'll be recording videos within minutes of opening the box.

Medical Grade

With its steel case and fluid resistant front panel, the USB200 is built to withstand the rigors of the OR, and it's so compact it will fit on any cart.

Features

- Captures and stores video and still digital images on USB flash drive.
- Color LCD monitor on front panel.
- Easy installation.
- Affordable, small, and portable.
- Powerful, friendly features.



Specifications

Dimensions

Width: 240 mm
 Height: 63 mm
 Depth: 200 mm
 Weight: 1.5 kg

Power

Voltage: 100 - 200 VCA
 Power Consumption: 15 W
 110/240 V3-prong Ac Jack

Temperature

Operation Temperature: -20°C to +40° C
 Storage Temperature: -40°C to +85° C

Format

Video format: MPEG2 (DVD-quality)
 Audio Format: MPEG1 L2, 2 Channels

Media

USB Flash Drives or External USB hard drive
 MediCapture brand USB flash drives (required)

Image format

Format: JPEG, TIFF, PNG, DICOM
 Image size: 1024 x 768, 800x 600, 640 x 480 Pixels

Audio

Inputs: Standard 3.5mm stereo phone jacks
 Outputs: Standard 3.5mm stereo phone jacks

Video

Input: 1xBNC, 1x S-Video (S-VHS, Y/C)
 Output: 1xBNC, 1x S-Video (S-VHS, Y/C)

Speakers

200 Hz-15kHz, 1W (for audio confirmation only)

LCD Monitor

2.5 in, full color (for Navigation Only)

Kit Contents

MediCap USB200
 USB flash drive
 video cables
 power cable
 User guide

