



Medivator Scope Buddy Endoscope Flushing Aid provides continuous, automated channel flushing, eliminating the risk of repetitive motion injuries. Simple one-touch operation quietly and quickly delivers fluids to endoscope channels. Compatible with all flexible endoscopes. Scope buddy can simultaneously flush multiple endoscopes shortening the time it would have taken to manually clean.

Features

- Replaces repetitive and labor-intensive syringing of channels.
- Hands-free operation reduces muscle fatigue and risk of injury.
- Simple push button operation.
- Continuous fluid delivery through all endoscope channels.
- Most endoscopes thoroughly flushed in 35 seconds.
- Simultaneously flushes multiple endoscope channels.
- Quickly and quietly delivers fluids through endoscope channels.
- Continuous, timed fluid delivery ensures consistent, repeatable endoscope flushing.
- Audible and visual alerts upon completion
- Compatible with all flexible endoscopes including ultrasound and those with elevator wire channels.



Specifications

Physical Specifications

9.5”H x 6”W x 8” D

Weight – 10 lbs

Electrical Requirements: 110 – 230 V, 50/60 Hz

Fluid Delivery Rate: 1 liter/minute

Flush Time: 0 – 99 minutes

System Pressure

Built-in safety bypass

Cycle Completion Indicator

Audible and visual alerts

Electrical Safety Certifications

UL Standard 61010-1

CANADA/USA C22.2 No. 61010-1

ETL and cETL listed

Unit Placement

Counter-top or wall-mounting
(bracket included)

Timer Control

Digital display with touchpad controls

